Selective Internal Radiation Therapy (SIRT)

PATIENT INFORMATION

Table of content

| | |
|--------------------------------------------------------------------------|----|
| Patient Information – English (English) | 3 |
| Informace o pacientovi – Čeština (Czech) | 6 |
| Patientinformation – Dansk (Danish) | 9 |
| Patienteninformation – Deutsch (German) | 12 |
| Información para el paciente – Español (Spanish) | 15 |
| Potilastiedot – Suomi (Finnish) | 18 |
| Informations à destination des patients – Français (French) | 21 |
| Informacije za pacijenta – Hrvatski (Croatian) | 24 |
| Betegtájékoztató – Magyar (Hungarian) | 27 |
| Informazioni per il paziente – Italiano (Italian) | 30 |
| პაციენტის ინფორმაცია – ქართული (Georgian) | 33 |
| Informacija pacientui – Lietuvių k. (Lithuanian) | 36 |
| Patiëntinformatie – Nederlands (Dutch) | 39 |
| Pasientinformasjon – Norsk (Norwegian) | 42 |
| Informacje dla pacjenta – język polski (Polish) | 45 |
| Folheto informativo do doente – Português (Portuguese) | 48 |
| Informații pentru pacienți – în limba română (Romanian) | 51 |
| Информация для пациента – русский (Russian) | 54 |
| Informácie pre pacientov – slovenčina (Slovak) | 57 |
| Informácie pre pacientov – Slovenčina (Slovenian) | 60 |
| Informacije za pacijenta – Srpski (Serbian) | 63 |
| Patientinformation – Svenska (Swedish) | 66 |
| Hasta Bilgileri – Türkçe (Turkish) | 69 |

Patient Information - English (English)

This document has been created for patients who are being considered to undergo or will undergo Selective Internal Radiation Therapy (SIRT) with QuiremScout™ and/or QuiremSpheres™. This document provides information about SIRT, and the products involved. Please contact your physician if you wish to receive additional details.


After the procedure your physician will give you an implant card. Please always carry that card with you and show it to any medical personnel who may be treating you:

What is SIRT?

SIRT is a therapy to treat liver cancer, where tiny radioactive beads (microspheres) are implanted into liver tumors. These liver tumors cannot be removed by surgery. SIRT is also known as Transarterial Radioembolization (TARE) or simply radioembolization.


What steps and products does SIRT include?


A typical SIRT therapy consists of two procedures. First comes the SIRT Work-up procedure, followed by the SIRT Treatment procedure.

How are QuiremScout™ and QuiremSpheres™ microspheres administered?

In both the 'SIRT Work-up' and 'SIRT Treatment' procedures, a small flexible tube known as a (micro) catheter is used. A radiologist places it in the body via the artery in the wrist or groin under local anesthesia. They then navigate the catheter to a pre-determined injection position in the liver. Microspheres are then administered into the blood stream leading to the tumor(s).


Liver arteries have a preferred flow towards liver tumors. Due to this, the microspheres flow directly towards the small blood vessels leading to the tumor. Due to the size of the microspheres, they become lodged into the small blood vessels. The implanted microspheres will stay there permanently and do not need to be removed.

How do QuiremScout™ and QuiremSpheres™ microspheres look like and what do they consist of?

The microspheres contained in QuiremScout™ and QuiremSpheres™ are identical in shape, size, and weight.

Microspheres are around 25 – 35 micrometers in diameter. This is equivalent to about a third of the width of a human hair.

A microsphere against a human hair.

Before being implanted, the microspheres are suspended in an administration liquid (2 ml).

Microspheres (pink) and administration liquid in glass vial.

QuiremScout™ and QuiremSpheres™ contain the radioactive element Holmium-166. The composition of both products is the following:

- Microspheres
Poly (L-lactic acid), Acetylacetone, Holmium chloride hexahydrate 99,9%
- Administration Liquid (water solution)
Di sodium hydrogen phosphate, Sodium di hydrogen phosphate, Kolliphor™ P188

How do QuiremScout™ and QuiremSpheres™ microspheres work?

The therapeutic effect of the microspheres is due to the radiation emitted by Holmium-166. The microspheres deliver:

- more than 90% of their radiation within the first 4 days, and
- near 100% within 8 days after the treatment procedure.

After implantation, microspheres can be viewed in the body via medical imaging (MRI, SPECT).

What should I remember before I undergo an MRI scan?

QiremScout™ and QiremSpheres™ have been assessed to be “MR conditional”. This means that you can safely undergo an MRI scan under certain conditions. By showing your implant card to the medical personnel who may be treating you, they will be able to determine these safe scanning conditions.

What radiation safety precautions should I take?

QiremScout™ and QiremSpheres™ are radioactive products. As a precaution, hospitals follow guidelines for the discharge of patients. These depend on local regulations and the dose you have received. The hospital will provide you with further information on these precautions when you leave the hospital.

Please inform your physician in case you are pregnant, breastfeeding or trying to conceive.

What side effects may I have?

Expected side effects are listed below. Contact your physician if you believe that you are experiencing side effects.

Expected side effects after ‘SIRT Work-up’ with QiremScout™

A common¹, moderate side effect is back pain. Less common², more severe side effects include abdominal³ and/or back pain.

Expected side effects after ‘SIRT Treatment’ with QiremSpheres™

Very common¹, mild to moderate side effects are:

| | | | |
|------------------------|-----------------------------------------------------|---------------------------------------|-------------------------------------------------------------------------------------------|
| Fatigue | Fever | Shortness of breath (dyspnea) | Abnormal accumulation of fluid in the abdominal cavity which may cause swelling (ascites) |
| Nausea and/or vomiting | Swollen limbs | Allergic reaction (to contrast agent) | |
| Dizziness | Irregular heartbeat (paroxysmal atrial tachycardia) | | |
| Back pain | | | |

Mild to severe abdominal pain can be expected.

Less common², more severe side effects are:

| | | | |
|---------|---------------|----------------------------------------------------------------|-----------------------|
| Nausea | Ascites | Bleeding in the esophagus ⁴ and/or digestive organs | Chest pain |
| Fatigue | Liver failure | | Ulcers in the stomach |

Laboratory blood values can change after SIRT with QiremSpheres™.

Possible serious side effects due to high radiation include inflammation of the :

| | | | | |
|-------------------------------|---------------------------|------------------------------------------------|-----------------------------|------------------------------|
| pancreas (acute pancreatitis) | stomach (acute gastritis) | gallbladder ⁵ (acute cholecystitis) | liver (radiation hepatitis) | lung (radiation pneumonitis) |
|-------------------------------|---------------------------|------------------------------------------------|-----------------------------|------------------------------|

and ulcers in the stomach and/or other digestive organs.

1. Common : Demonstrated in more than 1 out of 10 patients.
2. Less common : Demonstrated in 1 out of 10 to 1 out of 100 patients.
3. Abdominal pain : Pain felt between the chest and the groin, i.e. the region of the belly.
4. Esophagus : Tube between mouth and stomach.
5. Gallbladder : Small, pouch-like organ in the upper right part of belly which stores bile.

Manufacturer of QiremSpheres™ and QiremScout™:

Qirem Medical B.V.
Zutphenseweg 55, 7418 AH, Deventer
The Netherlands
qirem.com

Qirem Medical B.V. is a Terumo Company.

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Quirem Medical B.V.
Quirem Medical B.V. bir Terumo Şirketidir.