


Destination[™]

Guiding Sheath

GUIDING YOUR
PROCEDURES THROUGH
CONSISTENT RELIABILITY



Not all products are available for sale in all countries. This information is provided only in respect to markets where these products are approved or cleared. All products are not cleared or approved in the U.S.A. by the Food and Drug Administration. Please contact your Terumo local sales representative for more information.



MENU

ALLOWS **SMOOTH TRANSITION** AND EASY PENETRATION

SHAPES AND SIZES TO MEET YOUR PERIPHERAL
CHALLENGES

DESTINATION™ GUIDING SHEATHS



ALLOWS **SMOOTH TRANSITION** AND EASY PENETRATION

Atraumatic tip


- Minimizes potential for vessel damage

Smooth transitions

- Guidewire-to-dilator and dilator-to-sheath

Provides lesion access

- Exceptional trackability, even in the most challenging anatomy





SHAPES AND SIZES TO MEET YOUR PERIPHERAL CHALLENGES

Maintains shape and positioning

- Outstanding durability, even in the longest procedures

Enhances visualization for precise positioning


- Triple radiopacity (sheath, dilator, and gold coil marker)

Resists kinking

- Stainless steel coil for consistent reliability

Provides smooth movement

- Coil reinforced tubing and a PTFE inner layer that minimizes friction



- Cross Cut™ haemostatic valve (CCV)
- Tuohy-Borst valve (TBV)

CLICK TO
ZOOM

Designed to work with

NaviCross™ and **Glidewire Advantage™**

Support Catheter

Guide Wire

Data on file

Not all products are available for sale in all countries. This information is provided only in respect to markets where these products are approved or cleared. All products are not cleared or approved in the U.S.A. by the Food and Drug Administration. Please contact your Terumo local sales representative for more information.



SHAPES AND SIZES TO MEET YOUR PERIPHERAL CHALLENGES

Accommodates anatomical variations in vessel access

- Straight
- Hockey Stick
- Multi-Purpose
- Renal Double Curve (RDC)
- Left Internal Mammary Artery (LIMA)



- Cross Cut™ haemostatic valve (CCV)
- Tuohy-Borst valve (TBV)

CLICK TO
ZOOM

Designed to work with

NaviCross™ and **Glidewire Advantage™**



Support Catheter

Guide Wire

Data on file

Not all products are available for sale in all countries. This information is provided only in respect to markets where these products are approved or cleared. All products are not cleared or approved in the U.S.A. by the Food and Drug Administration. Please contact your Terumo local sales representative for more information.

ORDERING INFORMATION



Inner diameter Fr	Item reference	Length	Hydrophilic coating length	Distal curve	Valve
5FR 	GS*F5ST1C45	45 cm	5 cm	Straight	CCV
	GS*F5HS1C45	45 cm	5 cm	Hockey stick	CCV
	GS*F5MP1C45	45 cm	5 cm	Multipurpose	CCV
	GS*F5RDCC45	45 cm	5 cm	RDC	CCV
	GS*F5LIMC45	45 cm	5 cm	LIMA	CCV
	GS*K5ST1T90B	90 cm	15 cm	Straight	TBV
	GS*K5MP1T90B	90 cm	15 cm	Multipurpose	TBV
	GS*K5ST1C90B	90 cm	15 cm	Straight	CCV
	GS*K5MP1C90B	90 cm	15 cm	Multipurpose	CCV
6FR 	RSR01	45 cm	5 cm	Straight	CCV
	RSR02	45 cm	5 cm	Hockey stick	CCV
	RSR03	45 cm	5 cm	Multipurpose	CCV
	RSR13	45 cm	5 cm	RDC	CCV
	RSR14	45 cm	5 cm	LIMA	CCV
	RSP01	65 cm	35 cm	Straight	CCV
	RSC01	90 cm	15 cm	Straight	TBV
	RSC03	90 cm	15 cm	Multipurpose	TBV
	RSC05	90 cm	15 cm	Straight	CCV
	RSC07	90 cm	15 cm	Multipurpose	CCV

Packed 1 per box and includes dilator. All dilators are 0.038” wire compatible.

RDC: Renal Double Curve - LIMA: Left Internal Mammary Artery - CCV: Cross Cut Valve - TBV: Tuohy-Borst Valve
1.Data on file. Terumo Medical Corporation

Not all products are available for sale in all countries. This information is provided only in respect to markets where these products are approved or cleared. All products are not cleared or approved in the U.S.A. by the Food and Drug Administration. Please contact your Terumo local sales representative for more information.

ORDERING INFORMATION

Inner diameter Fr	Item reference	Length	Hydrophilic coating length	Distal curve	Valve
<div>7FR</div> <div></div>	RSR04	45 cm	5 cm	Straight	CCV
	54-74501	45 cm	35 cm	Straight	CCV
	RSR05	45 cm	5 cm	Hockey stick	CCV
	RSR06	45 cm	5 cm	Multipurpose	CCV
	RSR15	45 cm	5 cm	RDC	CCV
	RSR16	45 cm	5 cm	LIMA	CCV
	RSP02	65 cm	35 cm	Straight	CCV
	RSC02	90 cm	15 cm	Straight	TBV
	RSC04	90 cm	15 cm	Multipurpose	TBV
	RSC06	90 cm	15 cm	Straight	CCV
	RSC08	90 cm	15 cm	Multipurpose	CCV
<div>8FR</div> <div></div>	54-84501	45 cm	35 cm	Straight	CCV
	54-86501	65 cm	35 cm	Straight	CCV
	54-89001	90 cm	60 cm	Straight	CCV
	54-89006	90 cm	60 cm	Straight	TBV

Packed 1 per box and includes dilator. All dilators are 0.038” wire compatible.

RDC: Renal Double Curve - LIMA: Left Internal Mammary Artery - CCV: Cross Cut Valve - TBV: Tuohy-Borst Valve
1.Data on file. Terumo Medical Corporation

Not all products are available for sale in all countries. This information is provided only in respect to markets where these products are approved or cleared. All products are not cleared or approved in the U.S.A. by the Food and Drug Administration. Please contact your Terumo local sales representative for more information.



CONTACT

Terumo Corporation
+81 3 3374 8111

Terumo Europe NV
+32 16 38 12 11

**Terumo Interventional
Systems EMEA**
+33 1 47 16 09 30

EMEA Sales Offices

**Terumo Europe NV
Africa Business Division**
+32 16 38 13 08

**Terumo Europe NV
Benelux Sales Division**
Belgium:
0800 14468
The Netherlands:
0800 0231938

**Terumo Europe NV
Emerging Market Division**
+32 16 38 12 11

Terumo Deutschland GmbH
+49 6196 8023 0

**Terumo Deutschland GmbH
Österreich**
+43 2236 378020

**Terumo Deutschland GmbH
Switzerland**
+41 56 419 10 10

Terumo Europe España SL
+34 902 10 12 98

Terumo France S.A.S.
+33 1 30 96 13 00

Terumo Italia S.r.l.
+39 06 94 80 28 00

Terumo Russia LLC
+7 495 988 4740

Terumo Sweden AB
+46 3174 85 880

Terumo Middle East FZE
+971 4 292 0200

Terumo UK Ltd
+44 1276 480440

**Terumo BCT Tıbbi Cihazlar Dağıtım
ve Hizmetleri A.Ş.**
+902 16 645 92 00

Terumo Poland Sp. z o.o.
+48 22 120 16 00

 **Microvention Europe SARL**
30 bis, rue du Vieil Abreuvoir
78100 Saint-Germain-en-Laye France
+33(0)1 39 21 77 46



SHA

Mainta

- Out

Enhanc


- Trip

Resist

- Stai

Provid

- Coil
mini



• Cross Cut™ haemostatic valve (CCV)

• Tuohy-Borst valve (TBV)


Designed to work with
NaviCross™ and **Glidewire Advantage™**
Support Catheter Guide Wire



SHA

Accom

- Stra
- Hoc
- Mul
- Ren
- Left



Cross Cut™ haemostatic valve (CCV)

Tuohy-Borst valve (TBV)

Designed to work with
NaviCross™ and **Glidewire Advantage™**
Support Catheter Guide Wire